

# LC-MS Import and Public Viewer QA Fixture

**CONTROLLED QA TEST DOCUMENT**  
**NOT A REAL CERTIFICATE OR ANALYTICAL RESULT**

<b>Product name</b>	TEST Full Product Lambda 25MG
<b>Product code</b>	TEST-FULL-LAMBDA-25MG
<b>Batch / Lot</b>	QA-LOT-20260710-02
<b>Report type</b>	LC-MS
<b>Document purpose</b>	Bulk import, private storage, admin review, public visibility, and viewer QA only
<b>Scientific meaning</b>	No analytical measurements, purity claim, identity claim, certification, or approval is established

## QA Control Notes

- This PDF is intentionally non-sensitive and contains no real laboratory data.
- Filename, product code, batch/lot, and report type are deterministic for import matching.
- The imported report must default to Pending Review and Product Data Only unless an authorized administrator changes it.
- Public viewing must require explicit approval, explicit public enablement, a published product, and all existing file/security checks.
- A real LC-MS report should contain traceable instrument data, observed m/z, charge state, deconvoluted mass, and mass error. None are asserted here.

**Research Use Only - Controlled non-sensitive QA fixture**